

# Roadrunner

- ✓ Choose between 22", 24" or 26" wheels
- ✓ Hand-welded 7020 aluminium frame
- ✓ Designed and manufactured in Denmark
- ✓ Max. user weight 120 kg
- ✓ Weight from 12,5 kg

## What is Framerunning?

Framerunning is an innovative athletic discipline for disabled people with impaired balance. A Framerunner is a specially built 3-wheeled "running bike" without pedals, on which problems with balance is not an issue. Framerunning is a preferred sports discipline for people with cerebral palsy, arthritis or amputees.



### 1 POSITION

*Easy to adapt to the individual needs of the user thanks to the multi-adjustable saddle.*

### 2 HIGH-GRADE COMPENSATION

*It is easy to adjust the resistance of the handlebar when turning. The high-grade compensator. This increases the direction stability and helps keep the front wheel in a straight track.*

### 3 STABLE AND SAFE

*The steering angle on 55 degrees makes the Roadrunner more stable and more capable of handling involuntary/uncontrolled movements when moving fast.*



1



2



3



# Roadrunner for Framerunning

Strong and light framerunner that gives people with disabilities the freedom to run on their own.

The Roadrunner is a height-adjustable and stable running bike for young people and adults. The maximum user weight of the bike is 120 kg. Depending on the size of the wheels, the saddle height may vary between 70 and 95 cm from the ground.

The ultra lightweight frame weighs no more than 3,5 kg. If you include wheels, saddle, chest support and the handlebar the total weight of the Roadrunner is 12,5 kg.

The Roadrunner is onesize, but by changing the wheels into another wheel size, the bike can be adjusted to the preferred user height. That makes the Roadrunner perfect for the young user and also when he or she gets older. The Roadrunner is equipped with a disc brake on the front wheel as a standard. That makes it possible to change the wheel size without having to adjust other parts. The wheels are available in three sizes; 22", 24" or 26".

The Roadrunner is easy to adapt to the individual needs of the user, thanks to the multi-adjustable saddle and the chest support, which is adjustable in height, depth and angle as well as the handlebar, which is height- and depth adjustable. It is also possible to choose whether the brake should be forward or rear-pull, depending on your needs.



▲ White 24" Roadrunner  
with forward-pull brakes

MEASUREMENTS	22"	24"	26"
WHEEL SIZES	548 MM	600 MM	650 MM
HANDLEBAR   SADDLE	60 - 88 CM	60 - 88 CM	60 - 88 CM
HANDLEBAR   CHESTPLATE	20 - 60 CM	20 - 60 CM	20 - 60 CM
CHESTPLATE   SADDLE	10 - 45 CM	10 - 45 CM	10 - 45 CM
SADDLE HEIGHT	70 - 90 CM	72,5 - 92,5 CM	75 - 95 CM
CHESTPLATE HEIGHT	91 - 115 CM	93,5 - 117,5 CM	96 - 120 CM
HANDLEBAR WIDTH	61 CM	61 CM	61 CM
TOTAL WIDTH	82 CM	82 CM	82 CM
TOTAL LENGHT	185 CM	190 CM	195 CM
USER WEIGHT MAX	120 KG	120 KG	120 KG

