

Item Numbe	Quantity	Part Number	Part Name	Revisio	Comment
1	1	7350000250886	Spindel for 5 inch caste		
2	2	7350000250381	Lager 37 mm 25 mm		
3	1	7350000250244	M8X43 mm		
4	1	7350000250824	Sits medium		
5	1	7350000250893	Hjul axel 25 mm design		
6	1	7350000250541	Hjul 559 speed-runner		
7	1	7350000250213	Länkhjul 5 tum		
8	1	7350000250862	Fork gen2 Left		
9	1	7350000250909	Distance 4 inch caster		
10	2	7350000250947	Camberblock 25 mm 2,		
11	1	7350000250442	Hjul axel 25 mm design		
12	1	7350000250930	Ram Z-line medium låg		
13	1	7350000250367	Lager 28 12 mm		
14	1	7350000250374	Lager 28 15		
15	1	7350000250503	Drivring 25 tum		
16	1	7350000250923	Spännhållare Z-line		
17	1	7350000250916	Spänne Z-line		
18	1	7350000250428	M12		
19	1	7350000250954	Klamma fotstöd		
20	1		Fotbåge Z-line medium		
21	1	7350000250961	Hjul klamma Hope		